

Product	Size	REF
<b>Enamel</b>		
Enamel - E1, Extra-light	15g	615201
Enamel - E2, Light	15g	615202
Enamel - E3, Medium	15g	615203
Enamel - E4, Transparent	15g	615204
Enamel - E5, White	15g	615205
Enamel - E6, Dark	15g	615206
<b>Enamel Opal</b>		
Enamel Opal - EO1, Extra-light	15g	615211
Enamel Opal - EO2, Light	15g	615212
Enamel Opal - EO3, Medium	15g	615213
Enamel Opal - EO4, Transparent	15g	615214
Enamel Opal - EO5, HT	15g	615215
Enamel Opal - EO6, LT	15g	615216

<b>Enamel Effect</b>		
Enamel Effect - EE1, Sunrise	15g	615221
Enamel Effect - EE2, Violet	15g	615222
Enamel Effect - EE3, Sunset	15g	615223
Enamel Effect - EE4, Fog	15g	615224
Enamel Effect - EE5, Sky	15g	615225
Enamel Effect - EE6, Ivory	15g	615226
<b>Add-On Correction</b>		
Add-on Correction - C1, Light	15g	615401
Add-on Correction - C2, Medium	15g	615402
Add-on Correction - C3, Dark	15g	615403
Add-on Correction - C4, Transparent	15g	615404
<b>Add-On Gingiva</b>		
Add-on Gingiva - G1, Pink	15g	615411
Add-on Gingiva - G2, Reddish-pink	15g	615412
Add-on Gingiva - G3, Salmon	15g	615413
Add-on Gingiva - G4, Dark	15g	615414
Add-on Gingiva - G5, Mahogany	15g	615415

<b>Kit Configurations</b>		
Basic Kit (BL2, A1-A3.5, B1-B3, C1-C3, D2-D3)	EA	600601
CC Kit - Trial (A2, A3)	EA	600602
CC Kit - Aesthetic	EA	600603
CC Kit - Gingiva	EA	600604
Shade Series Kit	EA	600605
<b>Accessories</b>		
Shade Indicator - Basic 1	EA	D601002
Shade Indicator - Basic 2	EA	D601003
Shade Indicator - Aesthetic	EA	D601004
Shade Indicator - Gingiva	EA	D601005
Dentsply Sirona Modeling Liquid U	100mL	301126
Dentsply Sirona Modeling Liquid E	100mL	301127
Dentsply Sirona Modeling Liquid DU	50mL	650126
Ducera Liquid SD	50mL	5368272054
Ducera Liquid FORM	50mL	5368272034
Dentsply Sirona Stain and Glaze Liquid	15mL	601315

<b>Dentsply Sirona Universal Stain &amp; Overglaze</b>		
DS Universal Stain Kit (all items below)	EA	600700
DS Universal Body Stain - S0	5g	605520
DS Universal Body Stain - S1	5g	605521
DS Universal Body Stain - S2	5g	605522
DS Universal Body Stain - S3	5g	605523
DS Universal Body Stain - S4	5g	605524
DS Universal Incisal Stain - i1	5g	605531
DS Universal Incisal Stain - i2	5g	605532
DS Universal Stain - White	5g	605500
DS Universal Stain - Crème	5g	605501
DS Universal Stain - Sunset	5g	605502
DS Universal Stain - Copper	5g	605503
DS Universal Stain - Khaki	5g	605504
DS Universal Stain - Olive	5g	605505



New

Shade Series Dentins  
keyed to VITA®  
3D-Master shades

Product	3D-Master	REF
<b>Celtra Ceram, Shade Series Kit</b> Includes all items below		
Dentin 101	1M1	15g 615700
Dentin 102	1M2	15g 615701
Dentin 203	2L1.5	15g 615702
Dentin 204	2L2.5	15g 615703
Dentin 205	2M1	15g 615704
Dentin 206	2M2	15g 615705
Dentin 207	2M3	15g 615706
Dentin 208	2R1.5	15g 615707
Dentin 209	2R2.5	15g 615708
Dentin 310	3L1.5	15g 615709
Dentin 311	3L2.5	15g 615710
Dentin 312	3M1	15g 615711
Dentin 313	3M2	15g 615712
Dentin 314	3M3	15g 615713
Dentin 315	3R1.5	15g 615714
Dentin 316	3R2.5	15g 615715
Dentin 417	4L1.5	15g 615716
Dentin 418	4L2.5	15g 615717
Dentin 419	4M1	15g 615718
Dentin 420	4M2	15g 615719
Dentin 421	4M3	15g 615720
Dentin 422	4R1.5	15g 615721
Dentin 423	4R2.5	15g 615722
Dentin 524	5M1	15g 615723
Dentin 525	5M2	15g 615724
Dentin 526	5M3	15g 615725
Enamel - E6, Dark	n/a	15g 615206




Photos courtesy of Mario Pace, Kimmel Zahntechnik, GmbH, Koblenz, Germany

Dentsply Sirona | 570 West College Avenue | York, PA 17401 | 800-243-1942 | dentsplysirona.com/ceramics

DeguDent GmbH | Rodenbacher Chaussee 4 | 63457 Hanau-Wolfgang | Germany | +49 6181 59-50

© 2018 Dentsply Sirona Inc. All rights reserved. 22327/REV 2019-05



Celtra® Ceram  
All-Ceramic Veneering Porcelain

Versatility that delivers

Dentsply  
Sirona  
Lab

Dentsply  
Sirona  
Lab

\* VITA, VITA classical A1-D4, and 3D-MASTER are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co. KG. The product shades referred to in our comparative list are Dentsply Sirona's closest approximation to the referenced VITA 3D-Master shades.



# Versatility that delivers

Celtra Ceram introduces the versatility of a single veneering system designed exclusively for all-ceramic substructures including lithium silicates and zirconia oxides.

- + Versatile concept
- + Premium aesthetics
- + Consistent performance

Celtra Ceram is suitable for veneering all-ceramic frameworks. Coefficient of Thermal Expansion (CTE):  $9.0 \times 10^{-6} \text{ K}^{-1}$  (25-500°C)

- + Celtra Press, Zirconia-reinforced Lithium Silicate (ZLS) frameworks:  $9.7 \times 10^{-6} \text{ K}^{-1}$  (CTE 25-500°C)
- + Lithium Disilicate:  $10.0 - 10.5 \times 10^{-6} \text{ K}^{-1}$  (25-500°C)
- + Cercon® ht frameworks / Cercon Base:  $10.5 \times 10^{-6} \text{ K}^{-1}$  (25-500°C)
- + Cercon xt frameworks:  $10.1 \times 10^{-6} \text{ K}^{-1}$  (25-500°C)
- + Zirconia frameworks:  $10.1 - 11.0 \times 10^{-6} \text{ K}^{-1}$  (25-500°C)

## EXCEPTIONAL PERFORMANCE

For incredibly stable results



A proprietary feldspathic base formulation delivers an incredibly stable and robust result - every time, even after multiple firings

## SHADE CONFIDENCE

For lifelike, reproducible veneerings



Out-of-the-bottle shade accuracy: Fewer remakes, improved consistency



# Smart System

Celtra® Ceram masses reflect a smart system with basic and advanced powders offering versatility in the laboratory without the expense of comparable oversized systems.

### + Dentins

- 20 final shades (bleach/A-D)
- 26 Shade Series shades\*
- 6 Power Dentin (PD) shades
- 9 Dentin Effect (DE) colors
- 5 Dentin Gingiva (DG) shades
- 7 Opaqueous Dentin (OD) shades

### + Enamels


- 6 Enamel (E) shades
- 6 Enamel Opal (EO) shades
- 6 Enamel Effect (EE) colors

### + Add-ons

- 4 Correction (C) shades
- 5 Gingiva (G) shades

### Further Advantages

- Optimized particle size for smooth application and seamless buildup
- Correction/Add-on porcelain can be applied and fired with glaze/stain porcelain, or after glaze/stain firing



\*Keyed to VITA 3D-Master<sup>1</sup>

<sup>1</sup> VITA, VITA classical A1-D4, and 3D-MASTER are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co. KG.

## SINGLE VENEERING SYSTEM

For all-ceramic restorations



Compatible with both zirconia and lithium silicate frameworks

## UNIVERSAL FINISHING SYSTEM

For metal-based and all-ceramic restorations



Dentsply Sirona Universal Stains' broad compatibility eliminates multiple finishing systems in your lab - learn more at [dentsplysirona.com/ceramics](http://dentsplysirona.com/ceramics)

# Celtra® Ceram veneering porcelain



Celtra Ceram is a highly versatile, low-fusing veneering system designed exclusively for a broad range of all-ceramic substructures including lithium silicates and zirconia oxides. It is a smart system with expansive capabilities to meet the changing needs of laboratories worldwide.

### Versatile concept

- + Smart veneering concept suitable for both aesthetic and production environments utilizing core masses and multiple effect powders for all-ceramic restorations

### Premium aesthetics

- + Smart collection of core masses for high efficiency and yet a multitude of highly aesthetic effect porcelains satisfy the most creative minds in your dental laboratory

### Consistent performance

- + Shade accurate, stable and forgiveness achieved using a proprietary mineral source with seamless blending to all-ceramic substructures

## Ordering information

Product	Size	REF	Product	Size	REF
<b>Dentin</b>			Opaqueous Dentin - OD4	15g	615154
Dentin - BL1	15g	615130	Opaqueous Dentin - OD5	15g	615155
Dentin - BL2	15g	615131	Opaqueous Dentin - OD6	15g	615156
Dentin - BL3	15g	615132	<b>Dentin Effect</b>		
Dentin - BL4	15g	615133	Dentin Effect - DE1, Blue	15g	615161
Dentin - A1	15g	615134	Dentin Effect - DE2, Green-brown	15g	615162
Dentin - A2	15g	615135	Dentin Effect - DE3, Grey	15g	615163
Dentin - A3	15g	615136	Dentin Effect - DE4, Orange-brown	15g	615164
Dentin - A3.5	15g	615137	Dentin Effect - DE5, Pink	15g	615165
Dentin - A4	15g	615138	Dentin Effect - DE6, Violet	15g	615166
Dentin - B1	15g	615139	Dentin Effect - DE7, Yellow	15g	615167
Dentin - B2	15g	615140	Dentin Effect - DE8, White	15g	615168
Dentin - B3	15g	615141	Dentin Effect - DE9, Orange	15g	615169
Dentin - B4	15g	615142	<b>Dentin Gingiva</b>		
Dentin - C1	15g	615143	Dentin Gingiva - DG1, Pink	15g	615171
Dentin - C2	15g	615144	Dentin Gingiva - DG2, Reddish-pink	15g	615172
Dentin - C3	15g	615145	Dentin Gingiva - DG3, Salmon	15g	615173
Dentin - C4	15g	615146	Dentin Gingiva - DG4, Dark	15g	615174
Dentin - D2	15g	615147	Dentin Gingiva - DG5, Mahogany	15g	615175
Dentin - D3	15g	615148	<b>Power Dentin</b>		
Dentin - D4	15g	615149	Power Dentin - PD1	15g	615181
<b>Opaqueous Dentin</b>			Power Dentin - PD2	15g	615182
Opaqueous Dentin - ODO	15g	615150	Power Dentin - PD3	15g	615183
Opaqueous Dentin - OD1	15g	615151	Power Dentin - PD4	15g	615184
Opaqueous Dentin - OD2	15g	615152	Power Dentin - PD5	15g	615185
Opaqueous Dentin - OD3	15g	615153	Power Dentin - PD6	15g	615186