

+
**WE
KNOW
ENDO.**



MAILLEFER

EndoActivator®

An easy to use & safe solution
to activate your irrigation

Safe. Effective. Simple.

Contra-angled design for easy access to posterior teeth

Light-touch ON/OFF button

Ergonomic, cordless, and battery-operated handpiece



Quick connect/disconnect Activators for various procedures

3-speed sonic motor

Designed by Clifford J. Ruddle DDS, Robert H. Sharp DDS, Pierre Machtou DDS

One fact

Fluid activation, in well shaped canals, plays an important role in debridement and disinfection of the root canal system.

One solution

The EndoActivator® is designed to safely and vigorously energize intracanal irrigants during endodontic treatment.



Courtesy of Clifford J. Ruddle DDS

Multiple benefits

Safe

- Strong, flexible medical grade polymer tips
- Single patient use
- Uncoated and non-cutting tips

Effective

- Creates fluid hydrodynamics
- Improves debridement and the disruption of the smear layer and biofilm

Simple

- Very simple clinical technique
- Intuitive device
- Ideal in practice when portability is required

SYSTEM KIT INCLUDES	REF A 0911 000 000 00
1x EndoActivator® Handpiece/Driver (Battery Included)	REF A 0912 000 000 00
Activator Tips:	
25x Small (#15/.02)	X
25x Medium (#25/.04)	X
25x Large (#35/.04)	X
100 x Protective Barriers	REF A 0914 000 000 00
1x Directions for Use	X
1x Sequence Card	X
1x Instructional DVD	X
Possibility to order each part separately	

SPARE PARTS

	Small #15/.02 REF A 0913 022 015 00
	Medium #25/.04 REF A 0913 022 025 00
	Large #35/.04 REF A 0913 022 035 00

Technical modifications on our products are not subject to notifications. Photos of our instruments are not contractual. INL B EN ENAC BRC 000 02 2017 - updated 09/2019