

DECLARATION OF CONFORMITY



Distributor **W&H Sterilization S.r.l.**

Address Via Bolgara, 2
I-24060 Brusaporto (BG)
ITALY

Product **Multidem C17; Multidem C27;**

Model Code **19723111; 19723112**

Technical Specifications Demineralizer composed by Filter head , Filter cartridge, EB back flow preventer, Water block, tubing and fittings

Application field Laboratories and medicine; industrial.

Technical file TF 14

Compliances EN 1717:2000 (UNI EN 1717:2002) if equipped with:
> for fluid Category 2: EA back flow preventer (W&H optional cod. W320503X)
> for fluid Category 5: AB air gap (W&H optional cod. A812023X)

ISSUE DATE	APPROVED BY
28 March 2011	Simone Previcini General Management
W&H Sterilization S.r.l. - Via Bolgara, 2 - 24060 Brusaporto (BG) Italy	